



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton Micro Tissue Hand Homogenizer Tapered, 5 ml

€178.32
plus VAT &
Shipping

Product Images



Description

Tapered micro hand homogenizer, complete

For gradual comminution especially for soft tissues.

Coarse comminution takes place first in the upper part, fine in the lower part.

Additional Information

No.	WD-1365
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0956
EAN	4058072427726
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume	5 ml
Hand homogenizers and pestles features	can be clamped into stirrer pestle with ring/distance piece
Functionality	mechanical
Number of pistils	1
Material rod	Stainless steel
DM Mortar	13 mm
DM Rod	6.3 mm
DM Stamp	12.9 mm
Spare part	WD-1364 WD-1363
Shape stamp	conical. pointed
Length mortar	93 mm
Length pestle	220 mm
Material mortar	Borosilicate glass
Material stamp	Polytetrafluoroethylene (PTFE. Teflon)
Type homogenizer	Manual homogenizer. manual
Wall distance pestle 1	0.1 to 0.15 mm

