



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton Micro Tissue Hand Homogenizer Tapered, 15 ml

€205.33
plus VAT &
Shipping

Product Images



Description

Tapered micro hand homogenizer, complete

For gradual comminution especially for soft tissues.

Coarse comminution takes place first in the upper part, fine in the lower part.

Additional Information

No.	WD-1368
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0959
EAN	4058072427740
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume	15 ml
Hand homogenizers and pestles features	can be clamped into stirrer pestle with ring/distance piece
Number of pistils	1
DM Mortar	18 mm
DM Rod	6.3 mm
DM Stamp	17.9 mm
Spare part	WD-1367 WD-1366
Shape stamp	conical. pointed
Functionality	mechanical
Length mortar	114 mm
Length pestle	240 mm
Material mortar	Borosilicate glass
Material rod	Stainless steel
Material stamp	Polytetrafluoroethylene (PTFE. Teflon)
Type homogenizer	Manual homogenizer. manual
Wall distance pestle 1	0.1 to 0.15 mm

