



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton Micro Tissue Homogenizer w. Teflon® ring, 0.2 ml, 2 pcs./pack

€152.94
plus VAT &
Shipping

Product Images



Description

Tissue hand homogenizer

A Teflon® spacer straightens the mortar and reduces aerosol formation of the sample.

Chamber: 7 mm outer Ø, 33 mm length

PU: 2 pcs./pack

Additional Information

No.	WD-1301
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0906
EAN	4058072427412
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume	0.2 ml
Hand homogenizers and pestles features	Pestle with ring/distance piece
Functionality	mechanical
Number of pistils	1
Material rod	Borosilicate glass
DM Mortar	7 mm
DM Rod	6.3 mm
DM Stamp	6.9 mm
Shape stamp	conical. pointed
Length mortar	33 mm
Length pestle	115 mm
Material mortar	Borosilicate glass
Material stamp	Borosilicate glass
Type homogenizer	Manual homogenizer. manual
Wall distance pestle 1	0.1 to 0.15 mm

