



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Wheaton Micro Tissue Hand Homogenizer Dounce 1 ml, 2 Sets/Pack

**€328.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Dounce handheld homogenizers are ideal for enzyme studies.

They come with two pestles: one for coarser pre-grinding ("Loose") and one for complete homogenization ("Tight").

The spherical pestle ends protect cell nuclei in soft samples.

The mortar has a large reservoir and a pouring lip.

**Each set consists of 2 pestles + 1 mortar.**

**PU:** 2 Set/Pack

## Additional Information

No.	WD-1291
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0972
EAN	4058072427368
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume	1 ml
Functionality	mechanical
Number of pistils	2
Material rod	Borosilicate glass
DM Mortar	11 mm
DM Rod	6.3 mm
DM Stamp	10.9 mm
Shape stamp	around
Length mortar	48 mm
Length pestle	125 mm
Material mortar	Borosilicate glass
Material stamp	Borosilicate glass
Type homogenizer	Manual homogenizer. manual Dounce
Wall distance pestle 1	0.089 to 0.14 mm
Wall distance pestle 2	0.025 to 0.076 mm

