



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton Micro Tissue Hand Homogenizer Dounce 15 ml, 2 Sets/Pack

€412.74
plus VAT &
Shipping

Product Images



Description

Dounce handheld homogenizers are ideal for enzyme studies.

They come with two pestles: one for coarser pre-grinding ("Loose") and one for complete homogenization ("Tight").

The spherical pestle ends protect cell nuclei in soft samples.

The mortar has a large reservoir and a pouring lip.

Each set consists of 2 pestles + 1 mortar.

PU: 2 Set/Pack

Additional Information

No.	WD-1293
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0975
EAN	4058072427382
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume	15 ml
Functionality	mechanical
Number of pistils	2
Material rod	Borosilicate glass
DM Mortar	22 mm
DM Rod	6.3 mm
DM Stamp	21.9 mm
Shape stamp	around
Length mortar	94 mm
Length pestle	210 mm
Material mortar	Borosilicate glass
Material stamp	Borosilicate glass
Type homogenizer	Manual homogenizer. manual Dounce
Wall distance pestle 1	0.089 to 0.14 mm
Wall distance pestle 2	0.025 to 0.076 mm

