



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Wheaton Micro Tissue Hand Homogenizer Dounce 40 ml, 2 Sets/Pack

**€520.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Dounce handheld homogenizers are ideal for enzyme studies.

They come with two pestles: one for coarser pre-grinding ("Loose") and one for complete homogenization ("Tight").

The spherical pestle ends protect cell nuclei in soft samples.

The mortar has a large reservoir and a pouring lip.

**Each set consists of 2 pestles + 1 mortar.**

**PU:** 2 Set/Pack

## Additional Information

No.	WD-1294
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0976
EAN	4058072427399
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume	40 ml
Number of pistils	2
DM Mortar	32 mm
DM Rod	6.3 mm
DM Stamp	31.9 mm
Shape stamp	around
Functionality	mechanical
Length mortar	140 mm
Length pestle	285 mm
Material mortar	Borosilicate glass
Material rod	Borosilicate glass
Material stamp	Borosilicate glass
Type homogenizer	Manual homogenizer. manual Dounce
Wall distance pestle 1	0.089 to 0.14 mm
Wall distance pestle 2	0.025 to 0.076 mm

