



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton Potter-Elvehjem hand homogenizer with PVC coating, 10 ml, 2 pcs./pack

€328.00
plus VAT &
Shipping

Product Images



Description

Potter-Elvehjem homogenizers, complete.

In the event of glass breakage, the glass fragments are held in place by the safety coating, which also protects against mechanical impact, scratches or stains and provides particularly good grip.

Replacement pistils on request.

PU: 2 pcs./pack

Additional Information

No.	WD-1319
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0914
EAN	4058072427467
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume	10 ml
Hand homogenizers and pestles features	with splinter protection coating can be clamped in stirrer
Number of pistils	1
DM Mortar	16 mm
DM Rod	6.3 mm
DM Stamp	15.9 mm
Spare part	WD-1320
Shape stamp	around
Functionality	mechanical
Length mortar	74 mm
Length pestle	219 mm
Material mortar	Borosilicate glass
Material rod	Stainless steel
Material stamp	Polytetrafluoroethylene (PTFE. Teflon)
Type homogenizer	Manual homogenizer. manual Potter-Elvehjem
Wall distance pestle 1	0.1 to 0.15 mm

