



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton plunger for tissue hand homogenizer 9-0906, 2 pcs/pack

€183.24
plus VAT &
Shipping

Product Images



Description

Pestle with Teflon stamp

With stainless steel rod

Additional Information

| | |
|--|--|
| No. | WD-1357 |
| Manufacturer (Brand) | Wheaton |
| old neolab article no. | 9-0951 |
| EAN | 4058072427641 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Hand homogenizers and pestles features | can be clamped into stirrer |
| DM Rod | 6.3 mm |
| DM Stamp | 10.9 mm |
| Shape stamp | conical. pointed |
| Length pestle | 130 mm |
| Material rod | Stainless steel |
| Material stamp | Polytetrafluoroethylene (PTFE. Teflon) |
| for part no. | 9-0950 |
| for vessels with volume up to | 1 ml |

