



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® Headband Magnifier 2 x

€117.10
plus VAT &
Shipping

Product Images



Description

Comfort headband magnifier with binocular lens parts, enables binocular spatial vision. It is also suitable for spectacle wearers.

Additional Information

No.	2-9545
Manufacturer (Brand)	neoLab
EAN	4058072011086
Transport temperature	Room temperature
Color	Black
Magnifiers and glasses features	additional lenses/magnification levels available as accessories
Number of lenses / slices	2
Material lens	Polymethylmetacrylate (PMMA. PXM) (acrylic glass)
Working distance	250 mm
Material magnifier	Polyamide (PA. unspecified) glass fiber reinforced
Type magnifier	Headband magnifier
Magnification	2-fold

