



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Magnifier with frame, 2x

€143.10
plus VAT &
Shipping

Product Images



Description

Preset magnifier with lightweight frame, also suitable for non-eyeglass wearers.

Lens can be replaced at any time.

Delivery: frame with 1 lens part in box

Additional Information

No.	2-9540
Manufacturer (Brand)	neoLab
EAN	4058072011079
Transport temperature	Room temperature
Color	Black
Magnifiers and glasses features	additional lenses/magnification levels available as accessories
Number of lenses / slices	2
Material lens	Polymethylmetacrylate (PMMA. PXM) (acrylic glass)
Working distance	250 mm
Material magnifier	Aluminum anodized
Type magnifier	Front magnifier with bracket
Magnification	2-fold

