



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Magnifier lamp 1.75x biconvex, incl. illuminant 18 W + table clamp

€656.40
plus VAT &
Shipping

Product Images



Description

Daylight and distortion-free magnification.

Light source: energy-saving 18 W compact fluorescent lamp without disturbing heat radiation, incl. table clamp.

Lamp: approx. 8000 hrs. burning time

Lamp head: 3-D joint for individual adjustability in 3 axes

Illuminance: max. 6000 Lux at 30 cm distance

Lens: High quality, ground glass magnifier 171 x 108 mm, 3 diopters, magnification 175%, biconvex

Projection: max. 85 cm

Weight: 4.0 kg

Color : white

Additional Information

No.	3-1745
Manufacturer (Brand)	neoLab
EAN	4058072018214
Transport temperature	Room temperature
Color	White
Form	rectangular
Length	900 mm
Magnifiers and glasses features	with table clamp with illuminant with articulated stand with illumination
DM Lens	171 mm
Energy supply	Power grid
Power	18 W
Material lens	Soda-lime glass (normal glass) 1
Working distance	300 mm
Ground	biconvex
Illuminance	6000 Lux
Type illuminant	Fluorescent lamp
Spare part	3-1748
Material magnifier	Aluminum
Type magnifier	Illuminated magnifier with articulated stand
Magnification	1.75-fold

