



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® scale magnifier, magnification 10x

€138.30
plus VAT &
Shipping

Product Images



Description

Precision scale magnifier for various applications in laboratory, medicine and technology.

Distortion-free and razor-sharp image.

Easy-grip focusing ring. Achromatic glass lens, coated.

Supplied with standard scale (Fig. 1). Measuring range 20 mm, 1/10 graduation.

23 mm Ø , housing made of fiberglass-reinforced polyamide.

Supplied in rigid plastic box with space for other scales.

Additional Information

No.	2-9521
Manufacturer (Brand)	neoLab
EAN	4058072010959
Transport temperature	Room temperature
Color	White
Magnifiers and glasses features	with measuring scale with illumination
DM Lens	23 mm
Energy supply	Battery
Material lens	Polymethylmetacrylate (PMMA. PXM) (acrylic glass)
Ground	achromatic
Type illuminant	Incandescent lamp
Material magnifier	Polyamide (PA. unspecified) glass fiber reinforced
Type magnifier	Scale magnifier
Type measuring scale	Measuring range 20 mm. 1/10 graduation
Magnification	10-fold

