



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® stand magnifier round, 10 x
magnification (38 dpt), Ø lens: 35 mm,**

**€49.70
plus VAT &
Shipping**

Product Images



Description

Stand magnifier with a scratch-resistant aplanatic lens (\varnothing 35 mm) made of PXM[®] with a cera-tec[®] coating.

High magnification (38.0 dpt / 10x) for various applications. Optimal illumination from all sides is possible due to the transparent stand.

Additional Information

No.	2-2461
Manufacturer (Brand)	neoLab
EAN	4058072082338
Transport temperature	Room temperature
Color	transparent (crystal clear)
Form	around
DM Lens	35 mm
Material lens	Polymethylmetacrylate (PMMA. PXM) (acrylic glass)
Ground	aplanatic
Material magnifier	Polymethylmetacrylate (PMMA. PXM) (acrylic glass)
Type magnifier	Stand magnifier
Magnification	10-fold

