



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



WLD-TEC inoculation loops made of special stainless steel eyelet i.Ø 3 mm

€30.60
plus VAT &
Shipping

Product Images



Description

Inoculation eyelets made of special stainless steel wire Ø 0.6 mm PU = 10 pieces eyelet i.Ø 3 mm

Additional Information

No.	WT-0067
Manufacturer (Brand)	WLD-TEC
VGKL number	103076373
EAN	4066292446427
Gross weight	0.1
Transport temperature	Room temperature
Net weight	0.10000 kg
Material	Stainless steel
Country of origin	Germany
Region of origin	Thuringia
Packing width	0.14000 m
Packing height	0.05000 m
Packing depth	0.10000 m
Customs tariff number	84169000
DM Wire	3 mm

