



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® universal polyamide ties, 250 x 4.8 mm, 100 pcs./pack

€22.80
plus VAT &
Shipping

Product Images



Description

Disposable cable tie made of polyamide with grid webs.

High tensile strength and high chemical resistance, can also withstand thermal stress. For sealing and autoclaving contaminated sacks and bags.

Must be destroyed to open.

Color: transparent natural

Additional Information

No.	1-7112
Manufacturer (Brand)	neoLab
EAN	4058072045135
autoclavable at 121°C	No
Color	White
Material	Polyamide 6.6 (PA 6.6. Nylon)
Wide	4.8 mm
Length	250 mm
TBST MAX	85 °C
TBST MIN	-40 °C

