



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® silicone rubber loops, 30 pcs./pack

€94.30
plus VAT &
Shipping

Product Images



Description

Adjustable rubber bands made of silicone rubber, extremely resistant to chemicals and weathering. Elastic in temperature range from -46°C to 121°C, autoclavable. Each loop can be pulled wider or tighter by a small tension band as needed. Loops are 5 mm wide and can be extended to 10 cm Ø.

In bags of 30 pieces (10 each red, blue and yellow).

Additional Information

No.	6-2291
Manufacturer (Brand)	neoLab
EAN	4058072025663
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue red yellow
Material	Silicone
Wide	5 mm
Length	100 mm
TBST MAX	121 °C
TBST MIN	-46 °C

