



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex ejection cannula with sealing cap for Calibrex™ 525/530 dispensers, 2.5/5/10 ml, 90 mm

**€46.70**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Ejection cannulas for bottle-top dispensers

Can be replaced without tools; completely autoclavable.

Standard: length 90 mm for models up to 10 ml, 120 mm from 25 ml.

## Additional Information

---

No.	SR-0240
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7617
EAN	4058072367589
Transport temperature	Room temperature

