



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Ejection Cannula with Sealing Cap for Calibrex™ 525/530 Dispensers, 25 ml, 120 mm

**€46.70**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Ejection cannulas for bottle-top dispensers

Can be replaced without tools; completely autoclavable.

Standard: length 90 mm for models up to 10 ml, 120 mm from 25 ml.

## Additional Information

---

|                        |                  |
|------------------------|------------------|
| No.                    | SR-0241          |
| Manufacturer (Brand)   | Socorex          |
| old neoLab article no. | 7-7618           |
| EAN                    | 4058072367596    |
| Transport temperature  | Room temperature |

