



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Ejection Cannula with Cap for Calibrex™ 525/530 Dispensers, 50/100 ml, 120 mm

€46.70
plus VAT &
Shipping

Product Images



Description

Ejection cannulas for bottle-top dispensers

Can be replaced without tools; completely autoclavable.

Standard: length 90 mm for models up to 10 ml, 120 mm from 25 ml.

Additional Information

No.	SR-0242
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7619
EAN	4058072367602
Transport temperature	Room temperature

