



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Braun Petzold Sterican® cannula size 14 - G 23 x 1 ¼", 0.60 x 30 mm

€3.78
plus VAT &
Shipping

Product Images



Description

Braun Sterican® cannulas in various designs

- According to ISO 7864 and DIN 13097
- Thin-walled cannula
 - Cannula tube made of stainless chrome-nickel steel
 - With extremely smooth surface and finely dosed silicone coating
- Special faceted long/short ground joint
 - For low-pain punctures
- Transparent Luer-Lock plastic attachment
 - Made of polypropylene
 - Color coded according to ISO 6009
- Pack of 100 pieces

Additional Information

No.	TZ-1444
Manufacturer (Brand)	Braun Petzold
VGKL number	194210014
EAN	4058072386146
Transport temperature	Room temperature
diameter	0,6 mm
Material	Stainless steel
Length	30 mm
Manufacturer	Brown Petzold

