



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Braun Petzold Sterican® Cannula G 18 x 1½" SB, 1.20 x 40 mm

**€5.45**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Braun Sterican® cannulas in various designs

- According to ISO 7864 and DIN 13097
- Thin-walled cannula
  - Cannula tube made of stainless chrome-nickel steel
  - With extremely smooth surface and finely dosed silicone coating
- Special faceted long/short ground joint
  - For low-pain punctures
- Transparent Luer-Lock plastic attachment
  - Made of polypropylene
  - Color coded according to ISO 6009
- Pack of 100 pieces

## Additional Information

---

No.	TZ-1466
Manufacturer (Brand)	Braun Petzold
VGKL number	194211005
EAN	4058072980436
Transport temperature	Room temperature
diameter	1,2 mm
Material	Stainless steel
Length	40 mm
Manufacturer	Brown Petzold

