



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Drencher, straight, 0.9 x 25 mm

€26.00
plus VAT &
Shipping

Product Images



Description

Drenchers/Probes

Straight or curved drenchers with rounded extremities for targeted nasal or Mouth delivery in animals. Luer lock attachment for best fit on syringe. Made of stainless steel, autoclavable.

Additional Information

| | |
|------------------------|------------------|
| No. | SR-0290 |
| Manufacturer (Brand) | Socorex |
| old neolab article no. | 7-7691 |
| EAN | 4058072368074 |
| Transport temperature | Room temperature |
| diameter | 0,9 mm |
| Material | Stainless steel |
| Form | straight |
| Length | 25 mm |
| Cannulas properties | with Luer-Lock |

