



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Hypodermic Needles for Dosyst<sup>TM</sup> Self-Filling Syringes, 0.6 x 6 mm / 23G x 1/4"

**€16.74**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Needles and cannulas** in various designs.  
Luer lock made of stainless steel, autoclavable.

## Additional Information

---

No.	SR-0343
Manufacturer (Brand)	Socorex
old neolab article no.	7-7747
EAN	4058072368609
Transport temperature	Room temperature
diameter	0,6 mm
Material	Stainless steel
Form	straight
Length	6 mm
Cannulas properties	with Luer-Lock
Shape lace	oblique

