



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e



Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Needles for Dosys™ Self-Filling Syringes, 0.8 x 80 mm / 21G x 3 1/8"

**€34.74**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Needles and cannulas** in various designs.  
Luer lock made of stainless steel, autoclavable.

## Additional Information

---

No.	SR-0357
Manufacturer (Brand)	Socorex
old neolab article no.	7-7761
EAN	4058072368746
Transport temperature	Room temperature
diameter	0,8 mm
Material	Stainless steel
Form	straight
Length	80 mm
Cannulas properties	with Luer-Lock
Shape lace	oblique

