



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e



Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Needles for Dosyst<sup>TM</sup> Self-Filling Syringes, 1.8 x 10 mm / 15G x 3/8"

**€26.42**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Needles and cannulas** in various designs.  
Luer lock made of stainless steel, autoclavable.

## Additional Information

---

No.	SR-0379
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7783
EAN	4058072368968
Transport temperature	Room temperature
diameter	1,8 mm
Material	Stainless steel
Form	straight
Length	10 mm
Cannulas properties	with Luer-Lock
Shape lace	oblique

