



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Needles for DosystTM Self Filling Syringes, 2.0 x 30 mm / 14G x 1 1/8"

€21.42
plus VAT &
Shipping

Product Images



Description

Needles and cannulas in various designs.
Luer lock made of stainless steel, autoclavable.

Additional Information

No.	SR-0381
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7785
EAN	4058072368982
Transport temperature	Room temperature
diameter	2 mm
Material	Stainless steel
Form	straight
Length	30 mm
Cannulas properties	with Luer-Lock
Shape lace	oblique

