



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Needles for DosystTM Self-Filling Syringes, 1 pc, 1.2 x 10 mm / 18G x 3/8"

€2.08
plus VAT &
Shipping

Product Images



Description

Needles and cannulas in various designs.
Luer lock made of stainless steel, autoclavable.

Additional Information

No.	SR-0384
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7788
EAN	4058072369019
Transport temperature	Room temperature
diameter	1,2 mm
Material	Stainless steel
Form	straight
Length	10 mm
Cannulas properties	with Luer-Lock
Shape lace	oblique

