



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cannulas Luer-Lock connection, 1.0 x 100 mm, 3 pcs./pack

€17.30
plus VAT &
Shipping

Product Images



Description

Cannulas with Luer-Lock connection; female

Material: Neck brass chrome-plated, cannula stainless steel with bevel.

Reusable, corrosion resistant.

3 pieces/pack each

Additional Information

No.	2-3110
Manufacturer (Brand)	neolab
EAN	4058072000714
Transport temperature	Room temperature
diameter	1 mm
Material	Stainless steel
Length	100 mm
Cannulas properties	with Luer-Lock

