



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cannulas Luer-Lock connection, 1.0 x 150 mm, 3 pcs./pack

€21.90
plus VAT &
Shipping

Product Images



Description

Cannulas with Luer-Lock connection; female.

Material: Brass chrome-plated neck, stainless steel cannula with bevelled edge.

Reusable, corrosion resistant.

3 pieces/pack each

Additional Information

No.	2-3112
Manufacturer (Brand)	neoLab
EAN	4058072000738
Transport temperature	Room temperature
diameter	1 mm
Material	Stainless steel
Length	150 mm
Cannulas properties	with Luer-Lock

