



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cannulas Luer-Lock connection, 1.2 x 80 mm, 3 pcs./pack

€12.60
plus VAT &
Shipping

Product Images



Description

Cannulas with Luer-Lock connection; female.

Material: Brass chrome-plated neck, stainless steel cannula with bevelled edge.

Reusable, corrosion resistant.

3 pieces/pack each

Additional Information

No.	2-3115
Manufacturer (Brand)	neoLab
EAN	4058072000769
Transport temperature	Room temperature
diameter	1,2 mm
Material	Stainless steel
Length	80 mm
Cannulas properties	with Luer-Lock

