



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cannulas Luer-Lock connection, 1.5 x 100 mm, 3 pcs./pack

€19.80
plus VAT &
Shipping

Product Images



Description

Cannulas with Luer-Lock connection.; female.

Material: Brass chrome-plated neck, stainless steel cannula with bevel.

Reusable, corrosion resistant.

3 pieces/pack each

Additional Information

No.	2-3121
Manufacturer (Brand)	neoLab
EAN	4058072000813
Transport temperature	Room temperature
diameter	1,5 mm
Material	Stainless steel
Length	100 mm
Cannulas properties	with Luer-Lock

