



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Manifold, 4-channel

€116.58
plus VAT &
Shipping

Product Images



Description

Manifold - Multiple distributor

Practical distribution or aspiration of media in plates with 24 or 96 wells.

Tight fit on syringes with Luer and Luer Lock attachment.

Stainless steel with PTFE stops, autoclavable.

Additional Information

No.	SR-0167
Manufacturer (Brand)	Socorex
old neolab article no.	7-7438
EAN	4058072366865
Transport temperature	Room temperature
Material	Stainless steel
Form	T-shaped
Cannulas properties	with Luer-Lock

