



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Weighing spoon with shuttle function,  
black, small, PS, antistatic, 100 pcs/pack**

**€31.40  
plus VAT &  
Shipping**

## Product Images



## Description

---

The **neoLab weighing buckets with shuttle function** are a unique economical solution that combines the function of a bucket and a weighing shuttle in one product. The weighing process has been streamlined, now you can simply scan or remove the material in the bucket and place it directly on the scale.

This eliminates the need for a weigh boat or weighing paper, and saves you valuable lab space. Made of anti-static polystyrene that is water resistant and reduces the likelihood of materials sticking to the scoop

## Additional Information

---

No.	2-6641
Manufacturer (Brand)	neoLab
EAN	4058072155681
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
Material	Polystyrene (PS)
sterile	No
Wide	63 mm
Form	pear-shaped
Length	127 mm
Spoon properties	with handle

