



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® aluminum tweezers, 95 mm long,
straight, round tip**

**€11.10
plus VAT &
Shipping**

Product Images



Description

Completely seamless aluminum tweezers, anodized, in a patented and cured design.

Highest precision from Switzerland, anti-magnetic and stainless.

The lightweight among precision tweezers.

Total length: 95 mm

Additional Information

No.	3-1230
Manufacturer (Brand)	neoLab
EAN	4058072012601
Transport temperature	Room temperature
Material	Aluminum
Length	95 mm
Tweezers features	antimagnetic
Figure	1
Bend tweezers / scissors	straight. not curved
Surface/ Structure	smooth
Tip width	2 mm
Peak strength	0.3 mm
Type tweezers	Standard (open)
Type lace	around

