



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Mini tweezers, straight tip, tapered, 70 mm long

€19.60
plus VAT &
Shipping

Product Images



Description

Mini tweezers of the highest quality for very fine work.
Anti-magnetic and acid-proof.
Total length: 70 mm

Additional Information

No.	1-1051
Manufacturer (Brand)	neoLab
EAN	4058072031664
Material	Stainless steel
Length	70 mm
Tweezers features	with special chemical resistance antimagnetic
Figure	2
Bend tweezers / scissors	straight. not curved
Shape straight flat tweezers	lightly drawn in
Surface/ Structure	smooth polished
Tip width	0.3 mm
Peak strength	0.2 mm
Type tweezers	Flat tweezers
Type lace	pointed

