




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® forceps "grip" type 2A, flat tip, 120 mm long

€24.20
plus VAT &
Shipping

Product Images



Description

For more comfortable work!

The soft-touch coating ensures that the "grip" tweezers have a secure grip and are protected from the effects of heat or cold.

The soft special coating prevents hand fatigue.

Traditional Swiss craftsmanship, made of high-quality stainless steel, anti-magnetic and acid-resistant.

Each pair of tweezers is individually tested under a magnifying glass.

Soft touch coating is not solvent resistant.

Additional Information

| | |
|------------------------------|---|
| No. | 2-0203 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072047801 |
| Material | Stainless steel |
| Length | 120 mm |
| Tweezers features | with special chemical resistance antimagnetic |
| Figure | 4 |
| Bend tweezers / scissors | straight. not curved |
| Shape straight flat tweezers | lightly drawn in |
| Surface/ Structure | smooth |
| Tip width | 2 mm |
| Peak strength | 0.3 mm |
| Type tweezers | Flat tweezers |
| Type lace | around |

