



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® forceps "grip" type 3, straight, very fine tip, 120 mm long

**€21.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

For more comfortable work!

The soft-touch coating ensures that the "grip" tweezers have a secure grip and are protected from the effects of heat or cold.

The soft special coating prevents hand fatigue.

Traditional Swiss craftsmanship, made of high-quality stainless steel, anti-magnetic and acid-resistant.

Each pair of tweezers is individually tested under a magnifying glass.

Soft touch coating is not solvent resistant.

## Additional Information

---

No.	2-0204
Manufacturer (Brand)	neoLab
EAN	4058072047818
Material	Stainless steel
Length	120 mm
Tweezers features	with special chemical resistance antimagnetic
Figure	5
Bend tweezers / scissors	straight. not curved
Shape straight flat tweezers	not retracted ("triangular")
Surface/ Structure	smooth
Tip width	0.3 mm
Peak strength	0.2 mm
Type tweezers	Flat tweezers
Type lace	very pointed

