



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® forceps "grip" type 5, straight, short tip, 110 mm long

€24.90
plus VAT &
Shipping

Product Images



Description

For more comfortable work!

The soft-touch coating ensures that the "grip" tweezers have a secure grip and are protected from the effects of heat or cold.

The soft special coating prevents hand fatigue.

Traditional Swiss craftsmanship, made of high-quality stainless steel, anti-magnetic and acid-resistant.

Each pair of tweezers is individually tested under a magnifying glass.

Soft touch coating is not solvent resistant.

Additional Information

No.	2-0206
Manufacturer (Brand)	neoLab
EAN	4058072047832
Material	Stainless steel
Length	110 mm
Tweezers features	with special chemical resistance antimagnetic
Figure	7
Bend tweezers / scissors	straight. not curved
Shape straight flat tweezers	strongly indented
Surface/ Structure	smooth
Tip width	0.2 mm
Peak strength	0.2 mm
Type tweezers	Flat tweezers
Type lace	very pointed

