



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® forceps "ion" type 00, standard tip, straight, 115 mm

€28.80
plus VAT &
Shipping

Product Images



Description

For more hygiene in sensitive workplaces!

The special coating releases antimicrobial silver ions that neutralize bacteria, germs and other pathogens. No resistance is created in the process.

Traditional craftsmanship from Switzerland, made of high-quality stainless steel, anti-magnetic and acid-resistant. Each tweezer is individually tested under the magnifying glass.

Additional Information

No.	2-0211
Manufacturer (Brand)	neoLab
EAN	4058072047887
Material	Stainless steel
Length	115 mm
Tweezers features	with special chemical resistance antimagnetic
Figure	2
Bend tweezers / scissors	straight. not curved
Shape straight flat tweezers	not retracted ("triangular")
Surface/ Structure	smooth polished
Tip width	0.4 mm
Peak strength	0.5 mm
Type tweezers	Flat tweezers
Type lace	very pointed

