



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® precision tweezers with Xtel plastic tip, 125 mm long, very fine

**€101.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Stainless steel precision tweezers in Swiss quality, with tip made of Xtel plastic.  
Scratch-proof, silicone and chloride free.  
Reinforced with fiberglass, anti-magnetic. Clean room compatible.  
Heat resistant up to 350°C for short time.

## Additional Information

---

No.	3-1238
Manufacturer (Brand)	neolab
EAN	4058072012663
Transport temperature	Room temperature
Length	125 mm
Tweezers features	antimagnetic
Bend tweezers / scissors	straight. not curved
Surface/ Structure	smooth
Tip width	0.5 mm
Peak strength	0.4 mm
Type tweezers	Standard (open)
Type lace	very pointed

