



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Antistatic forceps, 3 mm, straight, with serration

€19.40
plus VAT &
Shipping

Product Images



Description

High quality tweezers made of antistatic vestamide with 20% glass fiber reinforcement.

Good stiffness, practically unbreakable, electrically conductive.

Excellent chemical resistance to acids and alkalis.

Temperature resistance: up to 180°C

Total length: 120 mm

Color: black

Additional Information

No.	2-1216
Manufacturer (Brand)	neoLab
EAN	4058072055844
Transport temperature	Room temperature
Material	Polyamide 12 (PA 12. Vestamid)
Length	120 mm
Tweezers features	with special chemical resistance with toothing anti-magnetic
Bend tweezers / scissors	straight. not curved
Surface/ Structure	smooth
Tip width	3.2 mm
Peak strength	0.4 mm
Type tweezers	Standard (open)
Type lace	square flattened

