



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Antistatic forceps, 1 mm, curved, without serration

**€19.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

High quality tweezers made of antistatic vestamide with 20% glass fiber reinforcement.

Good stiffness, practically unbreakable, electrically conductive.

Excellent chemical resistance to acids and alkalis.

Temperature resistance: up to 180°C

Total length: 120 mm

Color: black

## Additional Information

---

No.	2-1219
Manufacturer (Brand)	neoLab
EAN	4058072055875
Transport temperature	Room temperature
Material	Polyamide 12 (PA 12. Vestamid)
Length	120 mm
Tweezers features	with special chemical resistance antimagnetic
Bend tweezers / scissors	straight. not curved
Surface/ Structure	smooth
Tip width	1 mm
Peak strength	0.3 mm
Type tweezers	Standard (open)
Type lace	angled (strongly horizontal) pointed

