



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Ceramic tweezers straight, 125 mm, width 1 mm w. serration

€70.04
plus VAT &
Shipping

Product Images



Description

High-tech tweezers made of solid fine ceramics, extraordinarily wear resistant, temperature resistant up to 1000°C, electrically insulating, rustproof, antistatic, high chemical resistance.

Please handle these tweezers carefully, as ceramic can break.

Additional Information

No.	3-1201
Manufacturer (Brand)	neoLab
EAN	4058072012489
Transport temperature	Room temperature
Material	Ceramics
Length	125 mm
Tweezers features	with special chemical resistance with toothing anti-magnetic
Bend tweezers / scissors	straight. not curved
Surface/ Structure	smooth
Tip width	1 mm
Peak strength	0.06 mm
Type tweezers	Standard (open)
Type lace	pointed

