



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Millipore membrane filter tweezers with smooth tips

€24.73
plus VAT &
Shipping

Product Images



Description

Membrane filter tweezers for gentle handling of membrane filters, no risk of damaging the filter.
Highly polished stainless steel with beveled, non-serrated tips. Autoclavable.

Additional Information

No.	MP-1635
Manufacturer (Brand)	Millipore
old neoLab article no.	7-9050
EAN	4058072335649
Transport temperature	Room temperature
Material	Stainless steel
Tweezers features	antimagnetic
Bend tweezers / scissors	straight. not curved
Surface/ Structure	smooth
Type tweezers	Standard (open)
Type lace	around

