



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



RSG Membrane filter forceps, 105 mm, straight w. flattened ends

€7.85 As low as

€5.23

**plus VAT &
Shipping**

Product Images



Description

Tweezers made of stainless steel for gentle handling of membrane filters.
Rounded and flattened at the front.
Tweezers are 115 mm long, autoclavable and flame-retardant.

Order No. 1-6131

The standard filter tweezers must be pressed to hold the filter.

Additional Information

No.	RG-1283
Manufacturer (Brand)	RSG
old neolab article no.	1-6131
EAN	4058072338688
Transport temperature	Room temperature
Material	Stainless steel
Length	115 mm
Tweezers features	antimagnetic
Bend tweezers / scissors	straight. not curved
Surface/ Structure	smooth polished
Tip width	5.5 mm
Peak strength	0.5 mm
Type tweezers	Standard (open)
Type lace	square flattened angled (slightly vertical)

