



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® PBTP tweezers with fiberglass, 3 mm, straight, with serration

€16.90
plus VAT &
Shipping

Product Images



Description

High-quality tweezers made of polybutylene terephthalate (PBTP) with 20% fiberglass content.

Good rigidity and best chemical resistance.

High insulation values, temperature resistant up to 200°C.

Length: 120 mm

Additional Information

No.	3-1211
Manufacturer (Brand)	neoLab
EAN	4058072012533
Transport temperature	Room temperature
Material	Polybutylene terephthalate (PBTB) fiberglass
Length	120 mm
Tweezers features	with special chemical resistance with toothings anti-magnetic
Bend tweezers / scissors	straight. not curved
Surface/ Structure	smooth
Tip width	3.2 mm
Peak strength	0.4 mm
Type tweezers	Standard (open)
Type lace	square flattened

