



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® PBTP forceps with fiberglass, 1 mm, curved, without serration

**€16.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

High-quality tweezers made of polybutylene terephthalate (PBTP) with 20% fiberglass content.  
Good rigidity and best chemical resistance.  
High insulation values, temperature resistant up to 200°C.  
Length: 120 mm

## Additional Information

---

No.	3-1214
Manufacturer (Brand)	neolab
EAN	4058072012564
Transport temperature	Room temperature
Material	Polybutylene terephthalate (PBTB) fiberglass
Length	120 mm
Tweezers features	with special chemical resistance antimagnetic
Bend tweezers / scissors	straight. not curved
Surface/ Structure	smooth
Tip width	1 mm
Peak strength	0.3 mm
Type tweezers	Standard (open)
Type lace	angled (strongly horizontal) pointed

