



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® stainless steel tweezers, 120 mm long,
straight, very fine, anti-magnetic**

**€52.60
plus VAT &
Shipping**

Product Images



Description

Tweezers of highest precision made of stainless steel.

Extremely fine tips, tempered.

Shape: straight

Total length: 120 mm

Supplied in a plastic case.

Order No. 2-1038

Anti-magnetic, acid-resistant, with slide lock

Additional Information

No.	2-1038
Manufacturer (Brand)	neoLab
EAN	4058072049157
Material	Stainless steel
Length	120 mm
Tweezers features	with special chemical resistance with sliding closure anti-magnetic
Bend tweezers / scissors	straight. not curved
Surface/ Structure	smooth
Tip width	0.3 mm
Peak strength	0.2 mm
Type tweezers	Standard (open)
Type lace	very pointed

