



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Tweezers Set

€19.30
plus VAT &
Shipping

Product Images



Description

Practical set of tweezers for routine work. Consisting of 3 tweezers in a case.

- stainless
- Color 2 x silver, 1 x black
- straight and curved
- 90 and 95 mm long respectively

Tip type:

- 1 x pointed
- 2 x flat tip 3 and 4 mm

Additional Information

No.	3-0303
Manufacturer (Brand)	neolab
EAN	4058072011420
Transport temperature	Room temperature
Material	Stainless steel
Length	90 mm
Tweezers features	in set anatomically shaped
Number of parts in set	3
Bend tweezers / scissors	slightly vertically bent straight. not bent
Surface/ Structure	smooth
Type tweezers	Standard (open)
Type lace	pointed square flattened

