



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® tweezers electrophoretically coated

**€5.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Tweezers for all routine work. Slightly serrated tip

- stainless
- black, powder coated
- 95 mm long
- pointed

Weight: approx. 12 g

## Additional Information

---

No.	3-0300
Manufacturer (Brand)	neoLab
EAN	4058072011390
Transport temperature	Room temperature
Material	Stainless steel
Length	95 mm
Tweezers features	with toothings antimagnetic
Bend tweezers / scissors	straight. not curved
Surface/ Structure	smooth
Type tweezers	Standard (open)
Type lace	pointed

